

EXHIBIT 6

This file may not print or view at 100%.
Die lines and color breaks do not print.
“Holding Area” text, box and rules do not print.

Black

Cyan

Magenta

Yellow

PMS 186

die_cut

Approved according to
SHR_SOP 04.07.08_Att A

Category 2

The Signature will not be printed!


LAN schustn2 - Mar 09, 2023

IPM lyonm2 - Mar 09, 2023

PM/LAB sovacoee - Mar 09, 2023

GD/LAY cookd4 - Mar 09, 2023

PMS-8001598 CARTON
3.86X6.4 AUTO-LOCK



ACCUCHEK®
Guide

Roche

SUPPORT IS A CLICK OR CALL AWAY.

accu-chek.com
1-800-858-8072

Carefully read the User's Manual before use
For testing glucose in fresh capillary whole blood
For single patient use only
For in vitro diagnostic use
Store between -13–158 °F
Use between 10–90 % relative humidity

ACCUCHEK®
Guide

9:38am 5/11/22
106
mg/dL
7:38am 85 mg/dL

ACCUCHEK® Guide

**SIMPLY TAKE
JUST ONE STRIP**
Spill-resistant SmartPack® vial
lets you take one and spill none.

**SIMPLY USE
THE ENTIRE END**
Easy-fill test strips — place a
small drop of blood anywhere
along the end instead of on a
tiny spot.

**SIMPLY SEE,
DAY OR NIGHT**
Strip port light for testing, anytime.

Use of this device on multiple patients may lead to transmission of
Human Immunodeficiency Virus (HIV), Hepatitis C Virus (HCV),
Hepatitis B Virus (HBV), or other bloodborne pathogens.

Roche

ACCUCHEK® Guide

CONTAINS

- Accu-Chek Guide blood glucose meter
(Bluetooth-enabled, batteries included)
- Accu-Chek Softclix lancing device and
10 lancets
- Carrying case
- Easy-to-use instructions
Instrucciones en español incluidas

The Accu-Chek Guide test strips and the Accu-Chek Guide
control solutions are sold separately.

Roche

ACCUCHEK® Guide

Blood Glucose Monitoring System
For Single Patient Use Only

SIMPLY MORE

- Advanced accuracy¹
- Test with less pain²
- Backlit display
- Strip port light
- Strip ejector
- Connects with
the mySugr app

Roche

Simply Save
on future Accu-Chek products

Scan for savings options

Or visit accu-chek.com/save

mySugr®

Automatically log results, access testing
insights, and share data with your
healthcare professional on the go,
with the mySugr app.

Download on the
App Store

GET IT ON
Google Play

¹10/10 accuracy: 95% of measured blood glucose results shall fall within
±10 mg/dL of lab reference value for blood glucose concentrations
<100 mg/dL and within ±10% for blood glucose concentrations ≥100 mg/dL.
Data on file.

²8 out of 10 surveyed users of other lancing devices agreed that the Accu-Chek
Softclix lancing device was less painful than the lancing device they were using
at the time. Based on a 2020 Roche Diabetes Care, Inc. Lancing Market
Research Study.

Roche Diabetes Care, Inc.
9115 Hague Road
Indianapolis, IN 46256, USA
www.accu-chek.com

Meter assembled in India. Carrying case assembled in China.
Lancing device assembled in Germany.
Lancets assembled in Japan.

Roche

PROOF OF PURCHASE
Accu-Chek Guide Care Kit
Cat. No. 08453071001
NRC 65702-0729-10

ACCUCHEK®
Guide

Magenta stroked lines on the callout layer do not print.
Those lines indicate no print (keep out) areas.

8001598 10mm